

Product Datasheet

SFPQ Antibody [Out of stock] (orb1536000)

Description	SFPQ Antibody [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 578-707 of human SFPQ (NP_005057.1). MMIRQREMEEQMRRQREESYSRMGYMDPRERDMRMGGGGAMNMGDPYGSG GQKFPPLGGGGGIGYEANPGVPPATMSGSMMSGDMRTERFGQGGAGPVGGQ GPRGMGPGTPAGYGRGREEYEGPNKKPRF
Target	SFPQ
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 72kDa/76kDa, while the observed MW by Western blot was 115kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF, IHC, IHC-P (1:100), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.