

Product Datasheet

NTRK3 / TRKC Antibody (orb1536039)

Description	NTRK3 / TRKC Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human NTRK3 (NP_001007157.1). KTEINCRPPDDGNLFPLEGQDSGNSNGNASINITDISRNITSIHIENWRSLSHTLN AVDMELYTGLQKLTIKNSGLRSIQPRAFAKNPHLRYNLSSNRLTTLWSQLFQTLS LRELQLEQNFFNCSCDIRWMQLWQEQGEAKLNSQNLICINADGSQPLFRMNIS QCDLPEISVSHVNLTVREGDNAVITCNGSGSPLPDVDWIVTGLQSINTHQTNLN WTNVHAINLTLVNVTSEDNGFTLTCAENVVGMSNASVALT
Target	NTRK3 / TRKC
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	1.45 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 68kDa/91kDa/92kDa/93kDa/94kDa, while the observed MW by Western blot was 110kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.