

Product Datasheet

NRP1 / Neuropilin 1 Antibody (orb1536042)

Description	NRP1 / Neuropilin 1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant protein of human NRP1
Target	NRP1 / Neuropilin 1
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.