

Product Datasheet

RPSA / Laminin Receptor Antibody [Out of stock] (orb1536133)

Description	RPSA / Laminin Receptor Antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, IP, WB
Target	RPSA / Laminin Receptor
Concentration	0.5 mg/ml
Storage	Store at -20°C. Long term: aliquot and store at -20°C. Reconstituted: store at 2°C to 8°C for up to 3 weeks. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), IP (2 - 10 µg/mg lysate), WB (0.5 - 1 µg/ml)
Expiration Date	12 months from date of receipt.