

Product Datasheet

KRT20 / CK20 / Cytokeratin 20 Antibody (clone Ks20.8) (orb1536139)

Description	KRT20 / CK20 / Cytokeratin 20 Antibody (clone Ks20.8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Semi-purified human cytokeratin preparation
Target	KRT20 / CK20 / Cytokeratin 20
Preservatives	Filter sterilized and diluted in buffer, pH 7.5, rAlbumin, Sodium Azide
Storage	Store at 2-8°C. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Paraffin and frozen. Control tissue: Colon Carcinoma, Colon Mucosa, Bladder. Localization: Cytoplasmic.
Isotype	IgG
Clonality	Monoclonal
Dilution Range	IHC, IHC-P (1:50)
Expiration Date	12 months from date of receipt.