

Product Datasheet

OR7C1 Antibody (170-250 aa, Internal) (orb1536157)

Description	OR7C1 Antibody (170-250 aa, Internal)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P, WB
Immunogen	Synthesized peptide derived from human Olfactory receptor 7C1
Target	OR7C1
Preservatives	PBS, 0.02% Sodium Azide, 50% Glycerol, 0.5% rAlbumin
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA (1:20000), IF (1:200 - 1:1000), IHC, IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.