

Product Datasheet

HNRPA1 / HnRNP A1 Antibody (Internal) (orb1536167)

Description	HNRPA1 / HnRNP A1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Bovine, Gallus, Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IHC-P, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the central region of human hnRNP A1.
Target	HNRPA1 / HnRNP A1
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ICC (1:100 - 1:500), IF, IHC, IHC-P (1:100 - 1:200), IP (1:10 - 1:100), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.