

Product Datasheet

CD119 / IFNGR1 Antibody (C-Terminus) (orb1536168)

Description	CD119 / IFNGR1 Antibody (C-Terminus)
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-Terminus of human IFNGR1.
Target	CD119 / IFNGR1
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ICC (1:100 - 1:500), IF (1:100 - 1:500), IHC, IHC-P (1:100 - 1:200), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.