

Product Datasheet

PIK3CG / PI3K Gamma Antibody (Internal) (orb1536299)

Description	PIK3CG / PI3K Gamma Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the central region of human PI3K p110 gamma.
Target	PIK3CG / PI3K Gamma
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:100 - 1:200), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.