

Product Datasheet

MAPK7 / ERK5 Antibody (orb1536305)

Description	MAPK7 / ERK5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant protein corresponding to human ERK5.
Target	MAPK7 / ERK5
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:100 - 1:500), WB (1:1000 - 1:2000)
Expiration Date	12 months from date of receipt.