

Product Datasheet

JAG1 / Jagged 1 Antibody (aa33-250) (orb1536339)

Description	JAG1 / Jagged 1 Antibody (aa33-250)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant JAG1 (Gly33-Asp250) expressed in E. coli.
Target	JAG1 / Jagged 1
Preservatives	PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol. PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 0.05-2ug/ml; ICC 5-20ug/ml; IHCP 20 ug/ml; WB 1-5ug/ml;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.