

Product Datasheet

GRN / Granulin Antibody (aa21-120) (orb1536344)

Description	GRN / Granulin Antibody (aa21-120)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant GRN (Pro21-Ser120) expressed in E. coli.
Target	GRN / Granulin
Preservatives	PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol. PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 0.05-2ug/ml; ICC 5-20ug/ml; IHCP 20 ug/ml; WB 0.5-2ug/ml;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.