

Product Datasheet

TMPRSS2 / Epitheliasin Antibody (orb1536422)

Description	TMPRSS2 / Epitheliasin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-400 of human TMPRSS2 (NP_001128571.1). GLLWKFMGSKCSNSGIECDSSGTCINPSNWCDGVSHCPGGEDENRCVRLYGPN FILQVYSSQRKSWHPVCQDDWNENYGRAACRDMGYKNNFYSSQGIVDDSGST SFMKLNTSAGNVDIYKKLYHSDACSSKAVVSLRCIACGVNLNSSRQSRIVGGESA LPGAWPWQVSLHVQNVHVCGGSIITPEWIVTAAHCVEKPLNNPWHWTAFAGIL RQSFMFYGAGYQVEKVISHPNYDSKTKNNDIALMKLQKPLTFNDLVKP
Target	TMPRSS2 / Epitheliasin
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	1.729 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 53kDa/57kDa, while the observed MW by Western blot was 52kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:200), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.