

Product Datasheet

IL6ra antibody (orb153644)

Description	Goat polyclonal antibody to IL6ra
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Il6ra (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Mouse Spleen and Thymus lysates and in lysates of cell lines NIH3T3 and NSO gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	YLLQFQLRYRPVWSK
Clonality	Polyclonal
MW	50.4
Dilution Range	ELISA: 1:2000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.