

Product Datasheet

SKIV2L2 Antibody (orb1536440)

Description	SKIV2L2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human SKIV2L2 (NP_056175.3). MADAFGDELFSVFEGDSTTAAGTKKDKEKDKGKWKGPPGSADKAGKRFDGKL QSESTNNGKNKRVDVFEGTDEPIFGKKPRIEESITEDLSLADLMPRVKVQSVETV EGCTHEVALPAEEDYLPLKPRVGKAAKEYPFILDAFQREAIQCVDNNQSVLVS AH TSAGKTVCAEYAIALALREKQRVFTSPIKALSNQKYREMYEEFQDVGLMTGDVTI NPTASCLVMTTEILRSMLYRGSEVMREVAWVIFDEIHVMRDSERGVVWEETIILLP DNVHYVFLSATIPNARQFAEWICHLHKQPC
Target	SKIV2L2
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 117kDa, while the observed MW by Western blot was 123kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:200), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.