


Product Datasheet

LMO6 antibody (orb153648)

Description	Goat polyclonal to LMO6.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	LMO6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell line Daudi (calculated MW of 68.6kDa according to NP_006141.2). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	HPSRRRH YQCD
Clonality	Polyclonal
MW	68.6
Entrez	4007
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi lysate us...