

Product Datasheet

GRN / Granulin Antibody (aa234-419) (orb1536484)

Description	GRN / Granulin Antibody (aa234-419)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant GRN (Met234-Gly419) expressed in E. coli.
Target	GRN / Granulin
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 1:100-1:5000; ICC 1:50-500; IHCP 20 ug/ml; WB 1:50-400;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.