

Product Datasheet

CHAT Antibody (aa517-732) (orb1536505)

Description	CHAT Antibody (aa517-732)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant ChAT (Asp517-Pro732)
Target	CHAT
Preservatives	PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol. PBS, pH 7.4, 0.05% Proclin-300, 50% glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 0.05-2ug/mL; ICC 5-20ug/mL; IHCP 20 ug/ml; WB 1-5ug/mL;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.