

## Product Datasheet

### CD45RB Antibody (clone MEM-55) (orb1536517)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | CD45RB Antibody (clone MEM-55)                                   |
| <b>Species/Host</b>        | Mouse                                                            |
| <b>Reactivity</b>          | Human, Primate                                                   |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | FC, IHC, IHC-P, WB                                               |
| <b>Immunogen</b>           | Human thymocytes and T lymphocytes.                              |
| <b>Target</b>              | CD45RB                                                           |
| <b>Preservatives</b>       | PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide |
| <b>Concentration</b>       | 1 mg/ml                                                          |
| <b>Storage</b>             | Store at 2°C to 8°C. Do not freeze.                              |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG1                                                             |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Dilution Range</b>      | Flo (0.5 - 2 µg/ml), IHC, IHC-P (10 µg/ml), WB                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |