

Product Datasheet

HLA-DP/DR Antibody (clone MEM-136) (orb1536520)

Description	HLA-DP/DR Antibody (clone MEM-136)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, IP, WB
Immunogen	PHA-activated peripheral blood lymphocytes
Target	HLA-DP/DR
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo, IHC, IHC-P (10 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.