

Product Datasheet

SYNPO / Synaptopodin Antibody (aa97-300) (orb1536587)

Description	SYNPO / Synaptopodin Antibody (aa97-300)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant SYNPO (Asn97-Met300) expressed in E. coli.
Target	SYNPO / Synaptopodin
Preservatives	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 1:100- 1:5000; ICC 1:50-500; WB 1:50-400;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.