

Product Datasheet

BST2 Antibody (orb1536636)

Description	BST2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-161 of human BST2 (NP_004326.1). NSEACRDGLRAVMECRNVTHLLQQELTEAQKGFQDVEAQAATCNHTVMALMA SLDAEKAQGQKKVEELEGIEITTLNHLKLDASAEVERLRRENQVLSVRIADKKYYP SSQDSS
Target	BST2
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 18kDa/19kDa, while the observed MW by Western blot was 36kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:50 - 1:200), IHC, IHC-P (1:200), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.