

Product Datasheet

S100A12 Antibody (orb1536656)

Description	S100A12 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100A12 (NP_005612.1). MTKLEEHLEGIVNIFHQYSVRKGHFDLSKGELKQLLTKEANTIKNIKDKAVIDEI FQGLDANQDEQVDFQEFISLVAIALKAAHYHTHKE
Target	S100A12
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 10kDa, while the observed MW by Western blot was Refer to Figures.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:200), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.