

Product Datasheet

LAMP1 / CD107a Antibody (orb1536708)

Description	LAMP1 / CD107a Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human LAMP1 (NP_005552.3). CGKENTSDDPSLVIAFGGRGHTLTNFTNRNATRYSVQLMSFVYNLSDTHLFPNASSK EIKTVESITDIRADIDKKYRCVSGTQVHMNNVTVTLHDAITQAYLSNSSFSRGETR CEQDRPSPPTAPPAPPSPSPVPKSPSVDKYNVSGTNGTCLLASMGQLQLNLTYE RKDNNTTVTRLLNINPNKTSASGSCGAHLVTLELHS
Target	LAMP1 / CD107a
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 38kDa/44kDa, while the observed MW by Western blot was 120kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:50 - 1:200), IHC, IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.