

Product Datasheet

LAMB2 / Laminin Beta 2 Antibody (aa1579-1798) (orb1536775)

Description	LAMB2 / Laminin Beta 2 Antibody (aa1579-1798)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1579-1798 of human LAMB2 (NP_002283.3). DVRRAEQLLDARRARSWAEDEKQKAETVQAALAEAQAQGAIRGAVAD TRDTEQTLTYQVQERMAGAERALSSAGERARQLDALLEALKLKRAAGNSLAATAE ETAGSAQGRAQAEQLLRGPLGDQYQTVKALAERKAQGVLAQAQAEQLRDEA RDLLQAAQDKLQRLQELEGTYEENERALESKAAQLDGLARMRSVLQAINLQVQ IYNTCQ
Target	LAMB2 / Laminin Beta 2
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 195kDa, while the observed MW by Western blot was 196kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:200), WB (1:1000 - 1:2000)
Expiration Date	12 months from date of receipt.