

Product Datasheet

Integrin Alpha5 Beta1 Antibody (orb1536784)

Description	Integrin Alpha5 Beta1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Recombinant Human ITGA5 & ITGB1 Heterodimer protein
Target	Integrin Alpha5 Beta1
Preservatives	Sterile Solution PBS, 5% Trehalose. Sterile Solution PBS, 5% Trehalose
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	ELISA (1:1000 - 1:2000), IHC, IHC-P (1:50)
Expiration Date	12 months from date of receipt.