

Product Datasheet

SMG6 Antibody (orb1536796)

Description	SMG6 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1419 of human SMG6 (NP_060045.4). PDTMGKEMGSQEGTRLEDEEEDVVEDFEEDSEAEGSGGEDDIRELRAKKLALA RKIAEQRRQEKIQAVLEDHSQMRQMELEIRPLFLVPDTNGFIDHLASLARLLES RKYLIVVPLIVINELDGLAKGQETDHRAGGYARVVQEKARKSIEFLEQRFESRDSC LRALTSRGNELESIAFRSEDITGQLGNDDLLSCCLHYCKDKAKDFMPASKEEPI RLLREVLLTDDRNLRVKALTRNVPVRDIPAFLTWAQVG
Target	SMG6
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	2.35 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 37kDa/57kDa/160kDa, while the observed MW by Western blot was 210kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.