

Product Datasheet

IL17F Antibody (aa31-163) (orb1536826)

Description	IL17F Antibody (aa31-163)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant IL17F (Arg31-Gln163) expressed in E. coli.
Target	IL17F
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ELISA 1:100-1:5000; ICC 1:50-500; WB 1:50-400;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (1:250 - 1:2500)
Expiration Date	12 months from date of receipt.