

Product Datasheet

AP2A1 / AP2-Alpha Antibody (aa1-260) (orb1536840)

Description	AP2A1 / AP2-Alpha Antibody (aa1-260)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Recombinant human AP-2 complex subunit alpha-1 protein (aa1-260).
Target	AP2A1 / AP2-Alpha
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	0.31 mg/ml
Storage	Short term: -20°C; Long term: -80°C; Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IHC, IHC-P (10 µg/ml)
Expiration Date	12 months from date of receipt.