

Product Datasheet

MER / MERTK Antibody (aa840-999) (orb1536848)

Description	MER / MERTK Antibody (aa840-999)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Recombinant human Tyrosine-protein kinase Mer protein (aa840-999).
Target	MER / MERTK
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Short term: -20°C; Long term: -80°C; Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IHC, IHC-P (1:20 - 1:200, 10 µg/ml)
Expiration Date	12 months from date of receipt.