

Product Datasheet

ROBO2 Antibody (orb1536941)

Description	ROBO2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ROBO2 (NP_001122401.1). TALLKCKATGDPLPVISWLKEGFTFPGRDPRATIQEQGTLQIKNLRISDTGTCTCV ATSSSGETSWSAVLDVTESGATISKNYDLSDLPGPSKPQVTDVTKNSVTLWSQ PGTPGTLPASAYIEAFSQSVSNWQTVANHVKTTLTYVRGLRPNTIYLFMVRIN PQGLSDPSPMSDPVRTQDISPPAQGVDRQVQKELGDVLRVRLHNP
Target	ROBO2
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	1.43 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 130kDa/151kDa/153kDa, while the observed MW by Western blot was 175kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC-P (1:200), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.