

Product Datasheet

RAB5B Antibody (orb1537009)

Description	RAB5B Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant human Ras-related protein Rab-5B protein (1-215AA)
Target	RAB5B
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IHC-P (10 µg/ml), WB (1:1000 - 1:5000)
Expiration Date	12 months from date of receipt.