

Product Datasheet

ITGB7 / Integrin Beta 7 Antibody (orb1537042)

Description	ITGB7 / Integrin Beta 7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full length protein of human Integrin beta 7
Target	ITGB7 / Integrin Beta 7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.. Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Western blot analysis of Integrin beta 7 expression in HT29 (A), mouse lung (B) whole cell lysates.
Clonality	Polyclonal
Dilution Range	IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.