

Product Datasheet

CD79B / CD79 Beta Antibody (clone CB3-1) (orb1537064)

Description	CD79B / CD79 Beta Antibody (clone CB3-1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP, WB
Immunogen	Fraction of Ig-associated molecules isolated from Ramos B cells
Target	CD79B / CD79 Beta
Preservatives	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 2°C to 8°C. Do not freeze.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
Dilution Range	Flo (1 - 5 µg/ml), IHC-P (10 µg/ml), IP, WB
Expiration Date	12 months from date of receipt.