

Product Datasheet

NOTCH2NL Antibody (aa1-236) (orb1537155)

Description	NOTCH2NL Antibody (aa1-236)
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant NOTCH2NL (Met1-Asn236)
Target	NOTCH2NL
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
Concentration	1 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ICC 5-20ug/mL; WB 0.5-2ug/mL;
Clonality	Polyclonal
Dilution Range	IHC-P (10 µg/ml), WB
Expiration Date	12 months from date of receipt.