



Product Datasheet

ALB antibody (orb153719)


Description	Mouse monoclonal antibody to ALB. Serum albumin is the most abundant protein in plasma. It accounts ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human Albumin (NP_000468.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	69kDa
Uniprot ID	P02768
Entrez	213
Dilution Range	WB: 1:200-500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of rat brai...



Immunohistochemical staining of mouse br...