

Product Datasheet

CD95 / FAS Antibody (C-Terminus) (orb1537204)

Description	CD95 / FAS Antibody (C-Terminus)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-Terminus of human CD95.
Target	CD95 / FAS
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (1:100), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.