

Product Datasheet

**SYK Antibody (aa5-360, clone SYK-01)
(orb1537213)**

Description

SYK Antibody (aa5-360, clone SYK-01)

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant human SYK.
Target	SYK
Preservatives	Phosphate-buffered solution, pH 7.2, 0.09% sodium azide.. Phosphate-buffered solution, pH 7.2, 0.09% sodium azide.
Concentration	0.5 mg/ml
Storage	Store at 4°C. Avoid freezing.
Note	For research use only
Application notes	Further information: Immunohistochemistry: LS-B1652 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B1652 was determined to be 10 ug/ml.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (1:500)
Expiration Date	12 months from date of receipt.