

Product Datasheet

SDC1 / Syndecan 1 / CD138 Antibody (aa23-254) (orb1537312)

Description	SDC1 / Syndecan 1 / CD138 Antibody (aa23-254)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	Recombinant human Syndecan-1 protein (aa23-254).
Target	SDC1 / Syndecan 1 / CD138
Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.. PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Storage	Short term: -20°C; Long term: -80°C; Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IF, IHC-P (10 µg/ml)
Expiration Date	12 months from date of receipt.