

Product Datasheet

MAPT / Tau Antibody (phospho-T548) (orb1537332)

Description	MAPT / Tau Antibody (phospho-T548)
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-Terminus of human TAU.
Target	MAPT / Tau
Preservatives	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
Concentration	1 mg/ml
Storage	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (1:100), WB (1:500 - 1:1000)
Expiration Date	12 months from date of receipt.