

Product Datasheet

ENPEP / Aminopeptidase A Antibody (clone OTI1C4) (orb1537340)

Description	ENPEP / Aminopeptidase A Antibody (clone OTI1C4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Full length human recombinant protein of human ENPEP(NP_001968) produced in HEK293T cell.
Target	ENPEP / Aminopeptidase A
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol, 1% rAlbumin
Concentration	0.77 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo (1:100), IF (1:100), IHC-P (10 µg/ml), WB (1:500)
Expiration Date	12 months from date of receipt.