

Product Datasheet

LRP2 / Megalin Antibody (orb1537357)

Description

LRP2 / Megalin Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

Recombinant full length protein of human LRP2

Target

LRP2 / Megalin

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.. Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Concentration

1 mg/ml

Storage

Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.

Note

For research use only

Application notes

Further information: Western blot analysis of LRP2 expression in mouse lung (A) whole cell lysates. Immunohistochemical analysis of LRP2 staining in human breast cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Clonality

Polyclonal

Dilution Range

IHC-P (10 µg/ml), WB (1:500 - 1:2000)

Expiration Date

12 months from date of receipt.