

Product Datasheet

CD118 / LIF Receptor Alpha Antibody (orb1537376)

Description	CD118 / LIF Receptor Alpha Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Recombinant Mouse LIFR protein
Target	CD118 / LIF Receptor Alpha
Preservatives	0.2 µm filtered solution in PBS. 0.2 µm filtered solution in PBS
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA, IHC-P (1:100)
Expiration Date	12 months from date of receipt.