

Product Datasheet

CD45RO Antibody (clone 4KB5) (orb1537382)

Description	CD45RO Antibody (clone 4KB5)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Target	CD45RO
Preservatives	Filter sterilized and diluted in buffer, pH 7.5, rAlbumin, Sodium Azide
Storage	Store at 2-8°C. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Paraffin and frozen. Control tissue: Tonsil, Lymph Node. Localization: Membranous.
Isotype	IgG1,k
Clonality	Monoclonal
Dilution Range	IHC-P (1:50)
Expiration Date	12 months from date of receipt.