

Product Datasheet

CASP9 / Caspase 9 Antibody (Cleaved Asp353) (orb1537415)

Description

CASP9 / Caspase 9 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	A synthesized peptide derived from human Caspase 9 (Cleaved-Asp353).
Target	CASP9 / Caspase 9
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
Concentration	1 mg/ml
Storage	Upon receipt, store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: For western blots: incubate membrane with diluted antibody overnight in 5% w/v milk , 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ICC (1:100 - 1:500), IF (1:100 - 1:500), IHC-P (10 µg/ml), Peptide-ELISA (1:20000 - 1:40000), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.