

Product Datasheet

ACPP / PAP Antibody (orb1537455)

Description

ACPP / PAP Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full length protein of human TMPase
Target	ACPP / PAP
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.. Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Western blot analysis of TMPase expression in MCF7 (A) whole cell lysates. Immunohistochemical analysis of TMPase staining in human prostate cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.
Clonality	Polyclonal
Dilution Range	IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.