

Product Datasheet

ACE2 / ACE-2 Antibody (aa773-805) (orb1537467)

Description	ACE2 / ACE-2 Antibody (aa773-805)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Target	ACE2 / ACE-2
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Storage	Maintain refrigerated at 2°C to 8°C for up to 2 weeks. For long term storage aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Ig
Clonality	Polyclonal
Dilution Range	ELISA, IHC-P (1:100), WB (1:1000)
Expiration Date	12 months from date of receipt.