

Product Datasheet

ACPP / PAP Antibody (aa217-386) (orb1537479)

Description	ACPP / PAP Antibody (aa217-386)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant ACP3 (Ser217-Asp386) expressed in E.coli.
Target	ACPP / PAP
Preservatives	0.01 M PBS, pH 7.4, 0.05% ProClin 300, 50% Glycerol
Concentration	0.5 mg/ml
Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Note	For research use only
Application notes	Further information: ICC 5-20ug/mL; WB 0.5-5ug/mL;
Clonality	Polyclonal
Dilution Range	ICC (5 - 20 µg/ml), IHC-P (1:100), WB (0.5 - 5 µg/ml)
Expiration Date	12 months from date of receipt.