

Product Datasheet

NCALD / Neurocalcin Delta Antibody (orb1537544)

Description	NCALD / Neurocalcin Delta Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human NCALD (NP_114430.2). MGKQNSKLRPEVMQDLLESTDFTEHEIQEWYKGFLRDCPSGHLSEEFKKIYGN FFPYGDASKFAEHVFTFDANGDGTIDFREFIILSVTSRGKLEQKLKWAFSMYD LDNGYISKAEMLEIVQAIYKMVSSVMKMPDESTPEKRTEKIFRQMDTNRDGKL SLEEFIRGAKSDPSIVRLLQCDPSSAGQF
Target	NCALD / Neurocalcin Delta
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	9.32 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 22kDa, while the observed MW by Western blot was 18kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:50 - 1:100), IHC-P (5 µg/ml), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.