

Product Datasheet

SUV420H1 Antibody (Preservative Free) (orb1537557)

Description	SUV420H1 Antibody (Preservative Free)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthetic peptide corresponding to the center region of the Human SUV420H1
Target	SUV420H1
Preservatives	0.2 µm filtered solution in PBS. 0.2 µm filtered solution in PBS
Storage	Store at 2°C to 8°C for up to 1 month. Aliquot and store at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Applications should be user optimized.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC-P (1:100)
Expiration Date	12 months from date of receipt.